

Human C5a protein, His and Fc tag (active)

Cat. No. GTX00416-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


Complement Component 5a (C5a) is a component of the complement system which plays a key role in promoting migration and adherence of neutrophils and monocytes to vessel walls. C5a has been proven to be able to induce chemotactic migration of THP-1 cells. Therefore, chemotaxis assay used 24-well microchemotaxis system was undertaken to detect the chemotactic effect of C5a on the human monocytic cell line THP-1. Briefly, THP-1 cells were seeded into the upper chambers (100 µl cell suspension, 10⁶ cells/ml in RPMI-1640 with 0.5% FBS) and C5a (50 and 100 ng/ml diluted separately in serum free RPMI-1640) was added in lower chamber with a polycarbonate filter (8 µm pore size) used to separate the two compartments. After incubation at 37°C with 5% CO₂ for 2h, the filter was removed, then cells in low chamber were observed by inverted microscope at low magnification (×100) and the number of migrated cells were counted at high magnification (×400) randomly (five fields for each filter). C5a is able to induce migration of THP-1 cells.

Observed MW (kDa) 44 kDa.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris (pH8.0) and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris (pH8.0), 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** His-Tag and Fc-Tag; Thr678~Arg751 (NP_001726.2)**Expression System** HEK293 cells**Purity** > 95%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00416-pro Functional Assay Image


The chemotactic effect of GTX00416-pro Human C5a protein (active).

Upper: (A) THP-1 cells were seeded into the upper chambers and serum free RPMI-1640 without C5a was added in lower chamber, then cells in lower chamber were observed at low magnification ($\times 100$) after incubation for 3hrs.


(B) THP-1 cells were seeded into the upper chambers and 50 ng/ml C5a was added in lower chamber, then cells in lower chamber were observed at low magnification ($\times 100$) after incubation for 3hrs.

(C) THP-1 cells were seeded into the upper chambers and 100 ng/ml C5a was added in lower chamber, then cells in lower chamber were observed at low magnification ($\times 100$) after incubation for 3hrs.

Lower: The number of migrated cells were counted at high magnification ($\times 400$) randomly (five fields for each filter).


GTX00416-pro Image

SDS-PAGE analysis of GTX00416-pro Human C5a protein (active).


GTX00416-pro Image

WB analysis of GTX00416-pro Human C5a protein (active).



For full product information, images and publications, please visit our [website](http://www.genetex.com).