

Human Ovary (Normal) tissue slides

Cat. No. GTX03349

Applications	IHC-P, ISH
Species	Human

Package
5 slide

PRODUCT

Summary

Tissue was fixed in formalin immediately after excision, and embedded in paraffin. . At least one of the tissue slides from each lot was stained with H & E to ensure the quality. Tissue sections are 5 µm in thickness, and mounted on positively charged glass slides.

Tissue type

Normal

Applications

Application Note

Please bake slides at 60°C for 30 minutes before use.

Properties

Storage

Store at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).