

# Broad Spectrum Negative Control

**Cat. No. GTX03597**

**Applications** IHC-P, IHC-Fr

**Package**  
6 ml

## PRODUCT

### Summary

Negative controls are used to determine the specificity of an antibody. No staining should be observed when the primary antibody is omitted. Use GTX03597 Broad Spectrum Negative Control in place of the primary antibody with a section to evaluate nonspecific staining and allow better interpretation of staining results. This product is specifically designed to work with rabbit, mouse and goat antibodies. The Broad Spectrum Negative Control has been titrated to work with Polymer-based Secondary Systems in manual and automated IHC platforms.

## Properties

**Form** Liquid

**Storage** Store at room temperature.

**Purification** Purified IgGs

**Conjugation** Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).