

ImmunoHisto-Sealer

Cat. No. GTX03606

Application IHC-P

Package
30 ml

PRODUCT

Summary

ImmunoHisto-Sealer is a preservation solution that acts as a sealant to improve tissue morphology and stain over an extended period of time. When applied to the stained section, specimens can be permanently mounted without chromogens fading. Additional benefits of the ImmunoHisto-Sealer include: saves hands on time when mounting slides, avoids the use of toxic chemicals such as xylene and eliminates accumulation of hazardous waste. This product is designed for tissue specimens stained with peroxidase, alkaline phosphatase-based systems and various chromogens. After treating with GTX03606 ImmunoHisto-Sealer, water solution chromogens can be mounted in Permanent Mounting Media.

PROPERTIES

Form Liquid

Storage Store at room temperature.


Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03606 IHC-P Image

IHC-P analysis of human tonsil tissue using Cytokeratin H.M.W antibody HRP Blue with/without GTX03606 ImmunoHisto-Sealer.



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 19 Page 1 of 1