

10X Rapid Histo-Sealer

Cat. No. GTX03614

Applications	IHC-P	Package 10 ml


חם	$\overline{}$	\Box	Γ
PR		ועו	π

Summary

This product is designed for tissue specimens stained with peroxidase, alkaline phosphatase-based methods using a variety ofchromogens.

Properties	
Form	Liquid
Storage	Store at room temperature.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
11010	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03614 IHC-P Image

IHC-P analysis of human melanoma tissue using GTX12502 Melanoma antibody [PNL2] with GTX03614 10X Rapid Histo-Sealer.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 05 / 11 Page 1 of 1