

# ACE2 antibody [MSVA-919R] HistoMAX™

# Cat. No. GTX04425

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (1) Package 500 µl, 100 µl

## PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

Go to Normal Tissue Gallery

Go to Cancer Tissue Gallery

# **Applications**

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Note: Manual staining: Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note Highly recommended for IHC-P in human tissues.

Properties	
Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment (around aa705-805) of human ACE2
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 12 / 06 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



## GTX04425 IHC-P Image

IHC-P analysis of human testis tissue using GTX04425 ACE2 antibody [MSVA-919R] HistoMAX™. Leydig cells and spermatocytes of the testis show a particularly strong ACE2 immunostaining in the testis.



## GTX04425 IHC-P Image

IHC-P analysis of human renal cortex tissue using GTX04425 ACE2 antibody [MSVA-919R] HistoMAX™. A strong ACE2 immunostaining occurs in proximal tubuli of the kidney.



## GTX04425 IHC-P Image

IHC-P analysis of human clear cell renal cell carcinoma (ccRCC) tissue using GTX04425 ACE2 antibody [MSVA-919R] HistoMAX™.

Clear cell renal cell carcinoma strong strong membraous ACE2 positivity.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 12 / 06 Page 2 of 2