

Dual-tip Alcohol-Resistant Laboratory Marker

Cat. No. GTX05193

Package
12 pcs

PRODUCT

For laboratory use only.
Ink is resistant to 75% alcohol exposure.
Two marking pen widths: 1.0 mm & 0.4 mm.

Summary

Laboratory Use:
Features: Once dried, ink is resistant to 75% alcohol exposure without fading. Please note that decreased alcohol resistance may occur on some non-porous surfaces.
Applications: Suitable for marking test tubes, culture dishes/flasks, and related items in laboratories. Works on any surface.

Properties

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX05193 Image



For full product information, images and publications, please visit our [website](#).