

Trident 25X Tris Glycine Transfer Buffer

Cat. No. GTX16347

Applications

WB

References (1)

Package

1 liter

Applications

Application Note

For 1000ml 1xTransfer Buffer, add 40ml 25x, 200ml methanol to 760ml H₂O; mix well. The final 1x: Tris 12 mM, Glycine 96 mM, pH 8.3. Suggestion for Novex gels.

Properties

Form Liquid

Storage Store at room temperature. Please note this product has a 9 months shelf life.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).