

Trident Total Protein Extraction Kit

Cat. No. GTX16372
Application IP, ELISA, Immunoblotting

Package
20 test, 5 test

APPLICATION

Application Note


The use of protease inhibitors is not necessary prior to extraction. However if downstream application takes significant amounts of time or the protein extract will be stored for longer period of time, addition of protease inhibitors to lysis buffer is recommended.

PROPERTIES

Storage Store Native Lysis Buffer at 4°C and store Denaturing Lysis Buffer at room temperature.

Note For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Spin column-based protein extraction and cell fractionation technologies were developed by Invent Biotechnologies, Inc.

DATA IMAGES


GTX16372 Image


For full product information, images and publications, please visit our [website](https://www.genetex.com).