

# Trident Blue Prestained Protein Ladder

**Cat. No. GTX16376**

**Applications**      WB

**Package**  
500 µl

## Applications

### Application Note

3 µl or 5 µl per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively.

1.5~2.5 µl per well for general Western transferring.

Apply more for thicker (> 1.5 mm) or larger gel.

## Properties

**Form**                      Liquid


**Storage**                      Stable for up to 2 weeks at 25°C Stable for up to 3 months at 4°C. For long term storage, store at -20°C.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX16376 WB Image

Migration patterns of Trident Blue Prestained Protein Ladder in different electrophoresis conditions



For full product information, images and publications, please visit our [website](#).