

DAB Enhancer

Cat. No. GTX20675

Applications

IHC-P, IHC-Fr

Package

10 ml

PRODUCT

Summary

DAB Enhancer

Applications

Application Note

1. Add one drop (50µL) of DAB enhancer in 1mL of DAB incubation buffer. 2. Add DAB chromogen and Peroxidase Block to the above buffer. 3. Incubate tissue sections with the incubation medium.

Properties

Form Liquid

Buffer 10ml concentrated solution.

Preservative No preservative

Storage Store at room temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).