

Multiple Species Tissue Array - Small Intestine

Cat. No. GTX21529

Application	IHC-P	Package 5 slide

PRODUCT

The Multiple Species Tissue Array (MSTA) slides were designed to study protein expression patterns in different cells and tissues from multiple species. Tissue slices from three different species are mounted on each MSTA slide which can then be treated as a single histological slide for H&E staining, immunohistochemistry, or in situ hybridization. This format allows a rapid analysis of protein expression and localization across different species. MSTA slides can also be used to quickly determine the species reactivity of a given antibody.

Summary

GeneTex's contain 5 μ m sections of formalin-fixed, paraffin-embedded tissues (up to 1 cm x 1 cm) affixed to positively charged glass slides. The same type of tissues from 3 different species (Human, Mouse, and Rat) are arrayed on slides from left to right. Tissue diagnosis was performed by certified pathologists.

Tissue type Normal

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

	Suggested dilution	Recommended dilution
IHC-P Assay dependent	IHC-P	Assay dependent

Not tested in other applications.

PROPERTIES		
Storage	Cool dry at room temperature. 4°C for long term storage.	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.

Date 2024 / 05 / 16 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com