

# Human Matched Pair Tissue Array - Bladder Tumor

**Cat. No. GTX21548**

<b>Application</b>	IHC-P
<b>Species</b>	Human

**Package**  
5 slide

## PRODUCT

### Summary

The Matched Pair Paraffin Tissue (MPPT) slides are designed for identifying tumor-specific/metastasis genes or proteins. Slices from normal and malignant tissues are mounted on each MPPT slide which can then be treated as a single histological slide for H&E staining, immunohistochemistry, or in situ hybridization. This format allows a rapid analysis of protein expression and localization across normal and cancerous tissue.

GeneTex's tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to positively charged glass slides. The normal and carcinoma tissues are arrayed on slides from left to right. Tissue diagnosis was performed by certified pathologists.

<b>Tissue type</b>	Tissue pair
--------------------	-------------

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Not tested in other applications.

## PROPERTIES

<b>Storage</b>	Cool dry at room temperature. 4°C for long term storage.
----------------	--

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).