

FluoroGel mounting medium

Cat. No. GTX28214

Applications ICC/IF, IHC-P, IHC-Fr, IHC

References (15)

Package

20 ml

PRODUCT

Summary FluoroGel mounting medium - fluorescence

Applications

Application Note

Prior to use turn the bottle dropper tip to the open position. Close the dropper tip tightly when finished mounting to avoid evaporation and solidifying of mounting medium.

Properties

Form Liquid

Storage Keep as concentrated solution, aliquot and store at 4°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).