

TissueBond slide adhesive

Cat. No. GTX29909

Application IHC-P, IHC-Fr

Package
7 ml

PRODUCT

Summary TissueBond slide adhesive

APPLICATION

Application Note

Prepare of working solution:

- i) Working solution should be prepared just prior to use.
 - ii) Add the entire content of TissueBond in a Acetone compatible container.
 - iii) Add dry acetone to make a volume of 350 ml and stir well. This solution is enough for coating 500 slides.
- *Once the bottle has been opened, use it immediately. Prolonged storage is not recommended. Discard unused working solution.

Application Note:

- 1) Place slides in a metallic slide holder. Clean slides by immersing them in acetone for 1 minute.
- 2) Immerse slides in TissueBond working solution for 30 sec.
- 3) Remove and drain the slides. Rinse the slides thoroughly with deionized water three times. Change the water after every rinse. Avoid creating water bubbles by gentle agitation of slide holder before drying.
- 4) Air dry slides at room temperature. Once dried, slides can be used immediately or can be stored in a box at room temperature for future use.

PROPERTIES

Form Liquid

Storage Keep as concentrated solution, aliquot and store at 4°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).