

ImmunoHistoMount™

Cat. No. GTX30922

Application IHC

Package
30 ml

PRODUCT

Summary ImmunoHistoMount™

APPLICATION

Application Note

Apply 3-4 drops of mounting medium directly on top of the specimen and spread out evenly by tilting the slide back and forth or spread evenly with a 0.2 ml plastic pipette tip making sure the tissue is not touched. Excess medium can be removed by touching the edges of slide against a paper towel. Let stand at room temperature for about 1 to 2 hours. Or the slides can be heated at 37-40°C for 40-60 minutes. The slides are ready for visualization under a microscope.

PROPERTIES

Form Liquid

Storage This product is stored at room temp, however for long term storage 2-8°C is recommended. The pH of this product is more stable at 2-8°C. DO NOT FREEZE.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).