

# Organo/Limonene Mount™

**Cat. No. GTX30923**

**Applications** IHC

**Package**  
30 ml

## PRODUCT

**Summary** Organo/Limonene Mount™

**Applications**

## Application Note

Dehydrate tissue or cell smear slides with dehydrating agents e.g. 30 to 100% alcohol followed by xylene, toluene or any other dehydrating reagents. Add 2-3 drops of Organo mounting medium. Apply coverslip carefully, avoiding air bubbles. Visualize slide under a microscope. The slides can be dried by leaving at room temperature overnight or by heating at 37-40oC for several hours.

## Properties

**Form** Liquid

**Storage** Store at room temperature.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 02 / 01 Page 1 of 1