

# Immuno In Situ Mount™

**Cat. No. GTX30924**

## Applications

IHC

## Package

30 ml

## PRODUCT

### Summary

Immuno In Situ Mount™

### Applications

### Application Note

Apply 3-4 drops directly on the top of specimen and spread out evenly by tilting back and fourth or spread evenly with a 0.2 ml plastic pipette tip. Let stand at room temperature for about 60-90 minutes or heat slides at 37-40oC for 40-50 minutes.

## Properties

### Form

Liquid

### Storage

Store at room temperature.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).