

# Fluoroshield™ with DABCO

**Cat. No. GTX30927**

**Application** ICC/IF, IHC

**Package**  
20 ml

## PRODUCT

**Summary** Fluoroshield™ with DABCO

## APPLICATION

### Application Note

Apply 3-4 drops of mounting medium directly on top of the specimen and spread out evenly by tilting back and forth or spread evenly with a 0.2 ml plastic pipette tip making sure the tissue is not touched. Excess medium can be removed by touching the edges of the slide on a paper towel. Let stand at room temperature for about 5 minutes. Apply cover slip, carefully avoiding air bubbles.

## PROPERTIES

**Form** Liquid

**Storage** Store at 2-8°C, Protect from light, DO NOT FREEZE.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).