

Immuno Mount with DABCO™

Cat. No. GTX30930

Applications

IHC

Package

30 ml

PRODUCT

Summary

Immuno Mount with DABCO™

Applications

Application Note

Apply 2-3 drops of mounting medium and apply coverslip carefully without getting air bubbles.

Properties

Form

Liquid

Storage

Store at 2-8°C, Protect from light, DO NOT FREEZE.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).