

Immuno Mount with DABCO™

Cat. No. GTX30930

Application IHC

Package
30 ml

PRODUCT

Summary Immuno Mount with DABCO™

APPLICATION

Application Note

Apply 2-3 drops of mounting medium and apply coverslip carefully without getting air bubbles.

PROPERTIES

Form Liquid

Storage Store at 2-8°C, Protect from light, DO NOT FREEZE.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).