

# Immuno Automation Buffer™ (10X)

**Cat. No. GTX30931**

**Applications**

WB, ICC/IF, ELISA, IHC

References ( 5 )

Package

100 ml

**PRODUCT**

**Summary**

Immuno Automation Buffer™ 10X

**Applications**

**Application Note**

Dilute this 10X buffer 10 times (e.g. 90 ml of deionized or distilled water + 10 ml of this buffer), mix well, and store 1 X buffer solution at room temperature or 2-8°C for several weeks. Wash 3-5X with this buffer in between steps in all immunoassays protocol.

**Properties**

**Form**

Liquid

**Storage**

Room temperature. If stored at 2-8°C, some sodium chloride will be precipitated, if this occurs, bring the buffer to room temperature and mix well till all sodium chloride is dissolved.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).