

# Alkaline Phosphatase Enhancer

**Cat. No. GTX30944**

**Applications**

ICC/IF, IHC, ISH

**Package**

250 ml

**PRODUCT**

**Summary**

Alkaline Phosphatase Enhancer

**Applications**

**Application Note**

1. After the alkaline phosphatase reaction, wash slide with buffer, and distilled/ deionized water.
2. Add few drops of this ready to use Enhancing Reagent to cover the tissue/cell smear, for 1-2 minutes and discard. DO NOT WASH SLIDE.
3. Add few drops of ready to use alkaline phosphatase chromogen, follow instructions as outlined for these chromogens.

**Properties**

**Storage**

Store at room temperature.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).