

Triton X-100 (25%)

Cat. No. GTX30960

Applications

WB, ICC/IF, ELISA, IHC

References (1)

Package

50 ml

PRODUCT

Summary

Triton X-100 (25%)

Applications

Application Note

For diluting antibodies use 1:1000 in buffers (1 ml into 999 ml).

Washing buffer use 1:500 (1 ml into 499 ml).

The optimum dilution should be determined by the individual lab.

Properties

Form

Liquid

Buffer

This detergent is supplied as 25% in water.

Preservative

0.05% ProClin 300

Storage

Store at room temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).