

## Triton X-100 (25%)

## Cat. No. GTX30960

ApplicationWB, ICC/IF, ELISA, IHCPackage<br/>50 ml

PRODUCT

**Summary** Triton X-100 (25%)

## APPLICATION

## **Application Note**

For diluting antibodies use 1:1000 in buffers (1 ml into 999 ml).

Washing buffer use 1:500 (1 ml into 499 ml).

The optimum dilution should be determined by the individual lab.

PROPERTIES	
Form	Liquid
Buffer	This detergent is supplied as 25% in water.
Preservative	0.05% ProClin 300
Storage	Store at room temperature.
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2024 / 05 / 02 Page 1 of 1