

Trident Universal Protein Blocking Reagent (animal serum free)

Cat. No. GTX30963**Applications**

WB, ICC/IF, IHC-Fr, ELISA, IHC

References (11)

Package

100 ml

PRODUCT**Summary**

Protein Blocking Reagent (animal serum free)

Applications**Application Note**

Ready to use. Western blot: Immerse the filter in a sufficient amount of Blocking Reagent and incubate for 30 minutes at RT with gentle shaking. IHC/ICC: Incubate tissue sections/ cell smears in Blocking Reagent for 10 minutes at RT or 37°C. Note: Optimal incubation times should be determined by the researcher.

Properties**Form**

Liquid

Buffer

20mM Tris, 150mM NaCl, 0.02% Triton X-100

Preservative

0.05% ProClin 300

Storage

Store at 2-8°C.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).