

# Trident Universal Protein Blocking Reagent (animal serum free)

**Cat. No. GTX30963**
**Application** WB, ICC/IF, IHC-Fr, ELISA, IHC

Reference ( 12 )  
Package  
100 ml

## PRODUCT

**Summary** Protein Blocking Reagent (animal serum free)

## APPLICATION

### Application Note

Ready to use. Western blot: Immerse the filter in a sufficient amount of Blocking Reagent and incubate for 30 minutes at RT with gentle shaking. IHC/ICC: Incubate tissue sections/ cell smears in Blocking Reagent for 10 minutes at RT or 37°C. Note: Optimal incubation times should be determined by the researcher.

## PROPERTIES

**Form** Liquid

**Buffer** 20mM Tris 150mM pH7.4, Sodium Chloride, 0.02% Triton X-100

**Preservative** 0.05% ProClin 300

**Storage** Store at 2-8°C.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).