

Endogenous Alkaline phosphatase blocking kit

Cat. No. GTX30968

Application IHC

Package
15 ml

PRODUCT

Summary Endogenous Alkaline phosphatase blocking kit

APPLICATION

Application Note

After AP reaction GTX30968 is added to chromogen system as follows:

1. Chromogen buffer and substrates for AP are available as ready-to-use (one component), two components or three components. Mix according to manufacture instructions.
2. To one ml of this ready to use buffer, chromogen and substrate mixture, add two drop of GTX30968, mix well for few seconds.
3. Add this chromogen mixture to your specimen incubate according to manufactures instructions.

PROPERTIES

Storage Store at 2-8°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).