

Brij 35 (Polyoxyethylene Laural Ether)

Cat. No. GTx73318

Applications WB, IHC-P

Package
50 ml

PRODUCT

Summary Brij 35 (Polyoxyethylene Laural Ether) - 50ml

Applications

Application Note

For diluting antibodies use 1ml Brij 35 solution added to 400 ml of appropriate buffer. For washing buffer use 5 ml Brij 35 solution added to 95 ml of appropriate buffer. Optimal dilutions/concentrations should be determined by the individual lab.

Properties

Form Liquid**Buffer** This detergent is supplied as 15% in water**Preservative** 0.15 % Proclin 300**Storage** Store at room temperature.

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 01 / 31 Page 1 of 1