

# ImmunoStain blocker solution

**Cat. No. GTX73323**

**Applications**

IHC

**References** ( 4 )

**Package**

100 ml

## PRODUCT

### Summary

ImmunoStain Blocker solution ( A non-specific blocking reagent for immunocytochemical staining procedures)

### Applications

### Application Note

ImmunoStain blocker solution is ready to use and does not require any preparation.

Protocol Recommendations:

1. Deparaffinize slides and bring them to water.
2. If desired, block hydrogen peroxidase activity.
3. Wash tissue sections with wash buffer.
4. Apply enough drops of ImmunoStain blocker solution to cover tissue sections completely. Incubate at room temperature for 10 minutes.
5. Drain off excess of ImmunoStain blocker solution and apply primary antibody without washing slides(s).

## Properties

### Storage

Store at room temperature.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).