10X ImmunoStain wash buffer

Cat. No. GTX73344

Applications	IHC	References (1) Package 1 liter
PRODUCT		
Summary	10X ImmunoStain wash buffer (pH 7.5)	

Applications

Application Note

No pH adjustment or additives are required when using this formulation. Simply dilute the concentrated 10-fold buffer with distilled or deionized water. A total of 10 liters of working rinsing solution can be prepared with 1 liter of the concentrated reagent. This buffer is also the reagent of choice for rinsing tissue sections that are subjected to immunocytochemical staining using conventional horizontal processing techniques.

Properties		
Storage	Store at room temperature.	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
NOLE	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.