

Endo/Blocker

Cat. No. GTX73347

Applications	IHC	Package 50 ml

PRODUCT

This product is designed for use in peroxidase based immunohistochemical staining procedures done on paraffin tissue **Summary** sections to suppress nonspecific staining due to endogenous peroxidase activity.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC	Assay dependent

Note: Protocol Recommendations:

- 1. Incubate paraffin embedded issue sections for 5-10 minutes at room temperature prior to applying primary antibody.
- 2. Rinse tissue with wash buffer after incubation with peroxidase block and proceed to the next step.

Not tested in other applications.

Properties		
Storage	Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 02 Page 1 of 1