

# Endo/Blocker

**Cat. No. GTX73347**

**Application** IHC

**Package**  
50 ml

## PRODUCT

### Summary

This product is designed for use in peroxidase based immunohistochemical staining procedures done on paraffin tissue sections to suppress nonspecific staining due to endogenous peroxidase activity.

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

### Suggested dilution

### Recommended dilution

IHC

Assay dependent

### Note : Protocol Recommendations :

1. Incubate paraffin embedded issue sections for 5-10 minutes at room temperature prior to applying primary antibody.
2. Rinse tissue with wash buffer after incubation with peroxidase block and proceed to the next step.

Not tested in other applications.

## PROPERTIES

### Storage

Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).