



Fluorescence

Fluorescence
Immunostaining

Dylight conjugated anti-Mouse/Rabbit IgG antibodies

- High fluorescence intensity & photostability
- High signal-to-noise ratio
- Compatible with a broad pH range (pH 4-9)




Fluoroshield™ Mounting Medium with DAPI (GTX30920)

- Fortified with DAPI for DNA & RNA staining
- Preserve fluorescence during prolonged storage
- Prevent photobleaching
- Compatible with various fluorophore




Signal+ for Immunostaining (GTX49998)

- Enhances signals
- Reduces background
- Works with any systems
- Compatible with ABC or polymer complex
- Ready-to-use



Conventional buffer



Signal + Ab enhancer

PAP Pen (GTX22601)


- Keeps staining buffers localized
- Allows use of less reagents
- Insoluble in alcohol and acetone




Antibodies against Fluorescent Proteins



Anti-GFP antibody (GTX113617)



Anti-RFP antibody (GTX127897)



Anti-mCherry antibody (GTX128508)

Dylight Conjugated Secondary Antibodies

Citation Support

GTX213110-04	Rabbit IgG antibody (Dylight488)	Gt pAb	ICC/IF, IHC-Fr
GTX213110-05	Rabbit IgG antibody (Dylight594)	Gt pAb	ICC/IF, IHC-Fr
GTX76791	Rabbit IgG (H+L) antibody (DyLight 350)	Gt pAb	FACS, ICC/IF, IHC, IP, WB
GTX76793	Rabbit IgG (H+L) antibody (DyLight 680)	Gt pAb	FACS, ICC/IF, IHC, IP, WB
GTX76795	Rabbit IgG (H+L) antibody (DyLight 800)	Gt pAb	FACS, ICC/IF
GTX213111-04	Mouse IgG antibody (Dylight488)	Gt pAb	ICC/IF, IHC-Fr
GTX213111-05	Mouse IgG antibody (Dylight594)	Gt pAb	ICC/IF, IHC-Fr
GTX76786	Mouse IgG (H+L) antibody (DyLight 350)	Gt pAb	FACS, ICC/IF, IHC, IP, WB
GTX76787	Mouse IgG (H+L) antibody (DyLight 633)	Gt pAb	FACS, ICC/IF, IHC, IP, WB

Anti-Fluorescent Protein Antibodies

GTX113617	GFP antibody	Rb pAb	ELISA, ICC/IF, IHC, IHC-Fr, IHC-P, IP, WB, IHC-Wm
GTX628528	GFP antibody [GT859]	Ms mAb	ELISA, ICC/IF, IHC-Fr, IP, WB, IHC-Wm
GTX127897	RFP antibody	Rb pAb	ICC/IF, IP, WB
GTX628545	RFP antibody [GT1610]	Ms mAb	ICC/IF, IP, WB
GTX128508	mCherry antibody	Rb pAb	ICC/IF, IP, WB, IHC-Wm

Fluorescent Proteins

GTX65519	Enhanced Green Fluorescent Protein (EGFP) recombinant protein
GTX65495	Red Fluorescent Protein recombinant protein
GTX65497	Yellow Fluorescent Protein recombinant protein
GTX65494	Blue Fluorescence Protein (BFP) recombinant protein

Mounting Medium

GTX28214	FluoroGel mounting medium
GTX30920	Fluoroshield™ with DAPI
GTX30921	Fluoroshield™ with PI
GTX30927	Fluoroshield™ with DABCO

Counter Stain Chemical

GTX16206	DAPI Staining Solution
GTX14083	Propidium Iodide