

Human TNFRSF14 protein, His tag

Cat. No. GTX00145-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


TNFRSF14 (Tumor necrosis factor receptor superfamily member 14) belongs to the tumor necrosis factor receptor superfamily. TNFRSF14 functions in signal transduction pathways that activate inflammatory and inhibitory T-cell immune response. It binds herpes simplex virus (HSV) viral envelope glycoprotein D (gD), mediating its entry into cells. A binding ELISA assay was conducted to detect the association of TNFRSF14 with TNFα. Briefly, TNFRSF14 were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl TNFRSF14 were then transferred to TNFα-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-TNFRSF14 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of TNFRSF14 and TNFα was in a dose dependent manner.

Observed MW (kDa) 22 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Ser41~Ser209 (NP_003811.2)**Expression System** E. coli**Purity** > 95%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00145-pro Functional Assay Image

Functional ELISA analysis of GTX00145-pro Human TNFRSF14 protein which can bind immobilized TNF alpha protein.


GTX00145-pro Image

WB analysis of GTX00145-pro Human TNFRSF14 protein.


GTX00145-pro Image

SDS-PAGE analysis of GTX00145-pro Human TNFRSF14 protein.



For full product information, images and publications, please visit our [website](http://www.genetex.com).