

Human IL1F9 protein, His tag

Cat. No. GTX00148-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


IL36G (Interleukin 36 gamma), also known as IL-1F9 and IL-1H1 is a member of the IL-1 family. The receptor for IL36 gamma is reported to be a combination of IL-1 R6/IL-1 Rrp2 and IL-1 R3/IL-1 RAcP. Thus we have conducted a binding ELISA assay to detect the interaction of recombinant human IL36G with recombinant human IL1R1. Briefly, IL36G were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to IL1R1-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-IL1R1 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of IL36G with IL1R1 was in a dose dependent manner.

Observed MW (kDa) 19 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Met1~Ala162 (NP_001265497.1)**Expression System** E. coli**Purity** > 90%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx00148-pro Functional Assay Image

Functional ELISA analysis of GTx00148-pro Human IL1F9 protein which can bind immobilized IL1R1 protein.


GTx00148-pro Image

WB analysis of GTx00148-pro Human IL1F9 protein.


GTx00148-pro Image

SDS-PAGE analysis of GTx00148-pro Human IL1F9 protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).