

Human Vitamin D Receptor protein, His tag

Cat. No. GTX00163-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


Vitamin D Receptor (VDR) is the nuclear hormone receptor for vitamin D3. The receptor belongs to the family of trans-acting transcriptional regulatory factors and shows sequence similarity to the steroid and thyroid hormone receptors. It mediates the action of vitamin D3 by controlling the expression of hormone sensitive genes. Besides, E1A binding protein p300 (EP300) has been identified as an interactor of VDR, thus a binding ELISA assay was conducted to detect the interaction of recombinant human VDR and recombinant human EP300. Briefly, VDR were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to EP300-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-VDR pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of VDR and EP300 was in a dose dependent manner.

Observed MW (kDa) 20 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 100mM NaHCO₃ and 500mM NaCl (pH8.3) to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 100mM NaHCO₃ (pH8.3), 500mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Met272~Ser427 (NP_000367.1)**Expression System** E. coli**Purity** > 90%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00163-pro Functional Assay Image

Functional ELISA analysis of GTX00163-pro Human Vitamin D Receptor protein which can bind immobilized EP300 protein.


GTX00163-pro Image

SDS-PAGE analysis of GTX00163-pro Human Vitamin D Receptor protein.


GTX00163-pro Image

WB analysis of GTX00163-pro Human Vitamin D Receptor protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).