

Human Thymidine Kinase 1 protein, His and GST tag

Cat. No. GTX00185-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


Thymidine kinase 1 (TK1), is a human thymidine kinase. Thymidine kinase has been making a growing impact in the cancer research community. It has been found that elevated blood serum levels of TK-1 correlates with metastatic capabilities of the cancer and thereby can be used to detect malignant types of cancer, furthermore TK-1 has been found to show up in blood serum even before clinical symptoms even start to show. Besides, Glyceraldehyde 3-phosphate dehydrogenase (GAPDH) has been identified as an interactor of TK1, thus a binding ELISA assay was conducted to detect the interaction of recombinant human TK1 and recombinant human GAPDH. Briefly, TK1 were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to GAPDH-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-TK1 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of TK1 and GAPDH was in a dose dependent manner.

Observed MW (kDa) 55 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His and GST-Tag; confidential (NP_003249.3)**Expression System** E. coli**Purity** > 95%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx00185-pro Functional Assay Image

Functional ELISA analysis of GTx00185-pro Human Thymidine Kinase 1 protein which can bind immobilized GAPDH protein.


GTx00185-pro Image

WB analysis of GTx00185-pro Human Thymidine Kinase 1 protein.


GTx00185-pro Image

SDS-PAGE analysis of GTx00185-pro Human Thymidine Kinase 1 protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).