

Human Calpain 3 protein, His tag

Cat. No. GTX00194-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


Calpain 3 is a calcium-dependent cysteine protease mainly expressed in skeletal muscle. In humans, calpain 3 is encoded by the CAPN3 gene. This gene encodes a muscle-specific member of the calpain large subunit family that specifically binds to titin. Mutations in this gene are associated with limb-girdle muscular dystrophies type 2A. Alternate promoters and alternative splicing result in multiple transcript variants encoding different isoforms and some variants are ubiquitously expressed. Besides, Titin (TTN) has been identified as an interactor of CAPN3, thus a binding ELISA assay was conducted to detect the interaction of recombinant human CAPN3 and recombinant human TTN. Briefly, CAPN3 were diluted serially in PBS with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to TTN-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-CAPN3 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of CAPN3 and TTN was in a dose dependent manner.

Observed MW (kDa) 32 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Ile602~Ala821 (NP_000061.1)**Expression System** E. coli**Purity** > 95%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00194-pro Functional Assay Image

Functional ELISA analysis of GTX00194-pro Human Calpain 3 protein which can bind immobilized TTN protein.


GTX00194-pro Image

SDS-PAGE analysis of GTX00194-pro Human Calpain 3 protein.


GTX00194-pro Image

WB analysis of GTX00194-pro Human Calpain 3 protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).