

Human CYP3A4 protein, His tag

Cat. No. GTX00198-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note

Cytochrome P450 3A4 (CYP3A4) is an important enzyme in the body, mainly found in the liver and in the intestine. It oxidizes small foreign organic molecules (xenobiotics), such as toxins or drugs, so that they can be removed from the body. CYP3A4 is a member of the cytochrome P450 family of oxidizing enzymes and the CYP3A4 is the most common and the most versatile one in drug metabolism. Besides, Heat Shock 70kDa Protein 8 (HSPA8) has been identified as an interactor of CYP3A4, thus a binding ELISA assay was conducted to detect the interaction of recombinant human CYP3A4 and recombinant human HSPA8. Briefly, CYP3A4 were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to HSPA8-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-CYP3A4 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of CYP3A4 and HSPA8 was in a dose dependent manner.


Observed MW (kDa) 22 kDa.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Pro344~Asp497② (NP_001189784.1)**Expression System** E. coli**Purity** > 90%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX00198-pro Functional Assay Image

Functional ELISA analysis of GTX00198-pro Human CYP3A4 protein which can bind immobilized HSPA8 protein.



GTX00198-pro Image

SDS-PAGE analysis of GTX00198-pro Human CYP3A4 protein.



GTX00198-pro Image

WB analysis of GTX00198-pro Human CYP3A4 protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).