

Mouse IL2 Receptor alpha protein, His tag

Cat. No. GTX00298-pro

Applications Functional Assay**Species** Mouse**Package**

10 µg

Applications

Application Note


IL2RA (Interleukin-2 receptor alpha chain), also known as Tac antigen and as CD25, was initially identified as a 55kDa membrane glycoprotein which is capable of binding IL 2. Besides, mouse IL2 shares 75.1% AA sequence identity with rat IL2, suggesting the exist of cross-species activity. Thus we have conducted a binding ELISA assay to detect the interaction of recombinant rat IL2RA with recombinant mouse IL2. Briefly, IL2RA were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to IL2-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-IL2RA pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of IL2RA with IL2 was in a dose dependent manner.

Observed MW (kDa) 37 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal His-Tag; Glu22~Lys236 (NP_032393.3)**Expression System** E. coli**Purity** > 90%**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx00298-pro Functional Assay Image

Functional ELISA analysis of GTx00298-pro Mouse IL2 Receptor alpha protein which can bind immobilized IL2 protein.


GTx00298-pro Image

SDS-PAGE analysis of GTx00298-pro Mouse IL2 Receptor alpha protein.


GTx00298-pro Image

WB analysis of GTx00298-pro Mouse IL2 Receptor alpha protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).