

Mouse Perforin protein, His and GST tag

Cat. No. GTX00328-pro

Species Mouse

Package
10 µg

Applications

Observed MW (kDa) 65 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.

Preservative ProClin 300

Storage For short-term storage (1-2 weeks), store at 4°C. For long-term storage, store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.

Region/Sequence N-terminal His and GST-Tag; Val40~Lys355 (NP_035203.3)

Expression System E. coli


Purity > 90%

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES




GTX00328-pro Image

SDS-PAGE analysis of GTX00328-pro Mouse Perforin protein (active).



For full product information, images and publications, please visit our [website](https://www.genetex.com).



GTX00328-pro Image

WB analysis of GTX00328-pro Mouse Perforin protein (active).



For full product information, images and publications, please visit our [website](https://www.genetex.com).