

Rat IL15 protein (active)

Cat. No. GTX00385-pro

Applications Functional Assay

Species Rat

Package

10 µg

PRODUCT

Summary Recombinant rat IL-15 is a 13.5 kDa protein consisting of 115 amino acid residues.

Applications

Application Note

The ED₅₀ was determined by the dose-dependent stimulation of the proliferation of murine CTLL-2 cells is < 10 ng/ml, corresponding to a specific activity of > 1 x 10⁵ units/mg.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.

Region/Sequence MNWIDVRYDL EKIESLIQFI HDTTLYTDS DFHPSCKVTA MNCFLLELQV ILHEYSNMTL NETVRNVLYL ANSTLSSNKN VIESGCKECE
ELEERNFTEF LQSFHIVQM FINTS

Expression System E. coli

Purity > 98%

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).