

Mouse FGF9 protein

Cat. No. GTX00395-pro

Species Mouse

Package
10 µg

PRODUCT

Summary

Fibroblast Growth Factor-9 Mouse Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 205 amino acids and having a molecular mass of 23.308 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 10mM Tris HCl (pH8.0), 150mM Ammonium Sulfate.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The sequence of the first five N-terminal amino acids was determined and was found to be Pro-Leu-Gly-Glu-Val.

Expression System E. coli

Purity > 95% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).