

Mouse IL1 Beta protein

Cat. No. GTX00398-pro

Species Mouse

Package
10 µg

PRODUCT

Summary Interleukin-1 beta Mouse Recombinant produced in E.Coli is a non-glycosylated, Polypeptide chain containing 153 amino acids and having a molecular mass of 17.5 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O to 0.1-1.0mg/ml. Lyophilized from Sodium Citrate (pH4.0).

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MVPIRQLHYR LRDEQQKSLV LSDPYELKAL HLNQININQQ VIFSMFVQGEPSNDKIPVA LGLKGKNLYL SCVMKDGTPPT LQLESVDPKQ YPKKKMEKRFVFNKIEVKSK VEFESAEPN WYISTSQAEH KPVFLGNNSG QDIIDFTMES VSS.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).