

Rat IL1 Beta protein

Cat. No. GTX00410-pro

Species	Rat
----------------	-----

Package	10 µg
----------------	-------

PRODUCT

Summary Interleukin-1b Rat Recombinant produced in E.Coli is a non-glycosylated, Polypeptide chain containing 153 amino acids and having a molecular mass of 17.3 kDa.

Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	MVPIRQLHCRLRDEQKCLVSDPCELKALHLNGQNISQQWFSMSFVQGETSNDKIPVALGLKGLNLYLSCVMKDGTPTLQLESVDPKQYPKKKMEKRFVFNKIEVKTKVEFESAQFPNWWYSTSQAEHRPVFLGNSNGRDIIVDTMEPVSS.
Expression System	E. coli
Purity	> 95% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).