

Human Hsp90 alpha protein

Cat. No. GTX00441-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

Applications

Observed MW (kDa)	90 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	50mM Tris-HCl, 5mM Bme, 300mM NaCl, 10% Glycerol
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
---------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	Full-length protein without any tag; NP_001017963.2
-----------------	---

Expression System	E. coli
-------------------	---------


Purification	Affinity Purified
--------------	-------------------

Purity	> 90%
--------	-------

Conjugation	Unconjugated
-------------	--------------

Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	--

DATA IMAGES



GTX00441-pro Image

SDS-PAGE analysis of GTX00441-pro Human Hsp90 alpha protein.

For full product information, images and publications, please visit our [website](#).