

Human Hsp90 beta protein

Cat. No. GTX00442-pro

Applications Functional Assay**Species** Human**Package**

50 µg

Applications

Application Note


The protein has ATPase activity.

Observed MW (kDa) 90 kDa.

Properties

Form Liquid**Buffer** 20mM Tris, 175mM NaCl, 10% Glycerol, 1mM DTT, 0.1mM EDTA**Preservative** No preservatives**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full-length protein without any tag; NP_031381.2**Expression System** Baculovirus (Sf9 cells)**Purification** Affinity Purified**Purity** > 90%**Endotoxin** Endotoxin-free**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00442-pro Image

SDS-PAGE analysis of GTX00442-pro Human Hsp90 beta protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).