

Human Hsp70 protein, His tag (active)

Cat. No. GTX00448-pro

Applications Functional Assay**Species** Human**Package**

50 µg

Applications

Application Note

The protein tested positive for ATPase activity using a Malachite Green assay.


Observed MW (kDa) 70 kDa.

Properties

Form Liquid**Buffer** 20mM Tris buffer, 150mM NaCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full-length protein, His-Tag; NP_005336.3**Expression System** Baculovirus (Sf9 cells)**Purification** Affinity Purified**Purity** > 90%**Endotoxin** Endotoxin-free**Conjugation** Unconjugated**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX00448-pro Image

SDS-PAGE analysis of GTX00448-pro Human Hsp70 protein (active).



For full product information, images and publications, please visit our [website](https://www.genetex.com).