

Staphylococcus aureus Enterotoxin D antibody [B1913M]

Cat. No. GTX00453

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	ELISA, Sandwich ELISA
Reactivity	Staphylococcus aureus

Package
500 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX00454 ; Detection : GTX00453

Not tested in other applications.

Product Note No cross reaction to other types of Staphylococcus aureus enterotoxin A, B, C, and E.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.8 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Staphylococcus aureus Enterotoxin D
Purification	Protein A purified
Purity	> 95% pure (SDS-PAGE)
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 12 Page 1 of 1